

Nucleic Acid Isolation System

# QuickGene-Mini80

*One For Each Person*



# Personal nucleic acid isolation device, one for

QuickGene-Mini80 is a compact system requiring no centrifugation in the isolation process, giving less strain to samples and enabling rapid nucleic acid isolation. DNA/RNA can be easily isolated from various samples including whole blood/tissue/cells/plants/virus and others. The QuickGene-Mini80 shares a common isolation mechanism and kits with our semi-automated nucleic acid isolation system, QuickGene-810. QuickGene-Mini80 offers high performance comparable to QuickGene-810 at a reasonable price.


## Compact Design

- The small, lightweight QuickGene-Mini80 takes up minimal space on the lab bench and is easy to carry.
- QuickGene-Mini80 can be operated on the lab bench so, without moving around, the user can perform nucleic acid isolation.



## Revolutionary Porous Membrane

- The QuickGene-Mini80 uses patented porous membrane only 80µm thick.
- Because of the outstanding absorptive and desorptive performances of the membrane, high-purity nucleic acid can be easily obtained in high yield at low pressure.
- The ultra thin membrane enables nucleic acid isolation in shorter time than when compared to glass fiber membranes.




## High Purity, High Yield

- Highly reproducible QuickGene-Mini80 can stably isolate nucleic acid in high yield.
- The isolated DNA/RNA contains few impurities and can be used directly in PCR and RT-PCR.

### Length of isolated DNA

QuickGene DNA tissue kit S




M1: *A*Hind III digest  
M2: MidRange PFG Marker II  
M3: Mouse Genomic DNA(BD)

\*Balb/c mouse (female, 7week old)

### Isolation of plasmid DNA

QuickGene Plasmid kit SII



M : 1kb PLUS DNA Ladder


#### Yield and purity of plasmid DNA

|                 | Average* |
|-----------------|----------|
| DNA yield(µg)   | 23.0     |
| Purity(260/280) | 1.97     |
| Purity(260/230) | 2.28     |

\*Average of 4 samples

### Isolation of DNA from tissue

QuickGene DNA tissue kit S



M : 1kbp PLUS DNA Ladder


#### Yield and purity of genomic DNA

|                 | Tail  | Liver |
|-----------------|-------|-------|
| DNA yield(µg)   | 4.0   | 3.6   |
| Purity(260/280) | 1.916 | 1.925 |
| Purity(260/230) | 2.117 | 1.942 |

\*Balb/c mouse (female, 7week old)  
\*Without RNase treatment

### Isolation of RNA from tissue

QuickGene RNA tissue kit SII



M : 1kbp PLUS DNA Ladder

#### Yield and purity of total RNA

|                     | DNase treatment | Liver 32mg | Spleen 32mg | Heart 32mg |
|---------------------|-----------------|------------|-------------|------------|
| total RNA yield(µg) | +               | 155.4      | 71.8        | 21.2       |
|                     | -               | 159.3      | 74.3        | 21.6       |
| Purity(260/280)     | +               | 2.1        | 2.1         | 2.2        |
|                     | -               | 2.1        | 2.1         | 2.2        |
| Purity(260/230)     | +               | 2.3        | 2.1         | 2.0        |
|                     | -               | 2.2        | 2.1         | 1.8        |

\*Balb/c mouse (female, 7week old)

# each person

## Easy & Rapid Processing

- The operation is simple. Just set the sample and rotate the blue pressurizing Rotary Switch on both sides of the system.
- Centrifugation is not required for the isolation. The operation is simple, and saves valuable time.

### Easy and simple isolation operation




\*Commercially available microtubes are also usable.

### Workflow of isolation using DNA whole blood kit


DNA can be isolated from human whole blood in only 6 min.\*

\*Represents the actual time required to perform the operation according to our protocols.

#### Prepare lysate



#### Extraction process




#### pressurizing process

- ① Set holder into system
- ② Rotate pressurizing switch toward the front side to start pressurizing
- ③ Make sure that there is no residual liquid in the cartridge and return the pressurizing switch to original position
- ④ Pull out holder from system

### Isolation of RNA from cultured cells

QuickGene RNA cultured cell kit S



M :1kbp PLUS DNA Ladder


#### Yield and purity of total RNA

|                     | DNase treatment | Average* |
|---------------------|-----------------|----------|
| total RNA yield(µg) | +               | 12       |
|                     | -               | 11.1     |
| Purity(260/280)     | +               | 2.14     |
|                     | -               | 2.12     |
| Purity(260/230)     | +               | 2.17     |
|                     | -               | 2.25     |

\*Average of 4 samples

### Isolation of RNA from cells cultured in 6 cm dish

QuickGene RNA cultured cell HC kit S



M :1kbp PLUS DNA Ladder


#### Yield and purity of total RNA

| Cell type           | HL60   | HEK293*  | HeLa | COS-7 | NIH/3T3 |
|---------------------|--------|----------|------|-------|---------|
| Cell form           | Pellet | 6cm dish |      |       |         |
| total RNA yield(µg) | 45.6   | 110.5    | 46.9 | 43.0  | 34.4    |
| Purity(260/280)     | 1.99   | 1.86     | 2.19 | 2.07  | 1.93    |
| Purity(260/230)     | 2.01   | 1.98     | 1.98 | 1.96  | 1.91    |

\*HEK293 is the average of 4 samples.

### Isolation of RNA from blood cells

QuickGene RNA blood cell kit S



M :1kbp PLUS DNA Ladder

#### Yield and purity of total RNA

| Leukocyte(1x10 <sup>7</sup> ) | DNase treatment | Average* |
|-------------------------------|-----------------|----------|
| total RNA yield(µg)           | +               | 4.7      |
|                               | -               | 5.9      |
| Purity(260/280)               | +               | 2.30     |
|                               | -               | 2.10     |
| Purity(260/230)               | +               | 1.29     |
|                               | -               | 0.68     |

\*Average of 2 samples

## Easy isolation from various samples

Appropriate kit selectable depending on sample.

| Isolation example  | Purpose  | Corresponding kits                |
|--|--|-----------------------------------|
| <b>Whole blood, buffy coat</b><br>Human/ Cow, etc.: Whole blood (EDTA blood, heparin blood)  | genomic DNA isolation for genetic test, genotyping                               | <b>DNA whole blood kit S</b>      |
| <b>Tissue</b><br>Mouse/ Rat: Tail, Sperm, Liver<br><b>Materials, manufactured goods</b><br>Wheat/ Red bean/ Pleurotaceae/ Shimeji mushroom, etc. | DNA isolation for genome analysis of samples ranging from mouse/ rat to food     | <b>DNA tissue kit S</b>           |
| <b>Plasmid</b><br><i>E. coli</i>   | Miniprep for plasmid   | <b>Plasmid kit S II</b>           |
| <b>Tissue</b><br>Mouse/ Rat: Liver, Kidney, Brain, Spleen, Lung, Heart, Thymus   | total RNA isolation for expression analysis such as real-time PCR and RT-PCR     | <b>RNA tissue kit S II</b>        |
| <b>Cultured cell</b><br>Adherent/ Non-adherent (HeLa, HL60, etc.)<br><b>Tissue</b><br>Dicotyledonous plant (leaf)/ Tobacco (leaf), Plants        | total RNA isolation for expression analysis such as real-time PCR and RT-PCR     | <b>RNA cultured cell kit S</b>    |
| <b>Cultured cell</b><br>Cultured cell on 6cm/ 10cm dish  | High-yield total RNA isolation for northern blotting and microarray              | <b>RNA cultured cell HC kit S</b> |
| <b>Blood cell</b><br>Leukocytes separated by ammonium chloride<br>Erythrocytes of whole blood or Ficoll fraction                                 | High-purity total RNA isolation for expression profiling or various applications | <b>RNA blood cell kit S</b>       |

\*This chart shows isolation examples. Protocol may vary depending on samples and sample conditions.

## The specific set of isolation kits supports various samples

QuickGene isolation kits are optimized for the QuickGene system to isolate DNA and RNA in the shortest time and with the highest quality. Appropriate kit selectable depending on sample.

| Isolation kits   | Isolation example   |
|--|---|
| QuickGene DNA whole blood kit S [ for 96 samples ] Reference code DB-S       | ca.5 $\mu$ g / Whole blood 200 $\mu$ l                      |
| QuickGene DNA tissue kit S [ for 96 samples ] Reference code DT-S            | ca.4 $\mu$ g / 5mg Balb/c Mouse tail                        |
| QuickGene Plasmid kit S II [ for 96 samples ] Reference code PL-S2           | ca.12.5 $\mu$ g / 1ml culture pBlueScript II / DH5 $\alpha$ |
| QuickGene RNA tissue kit S II [ for 96 samples ] Reference code RT-S2        | ca.100 $\mu$ g / 30mg Mouse liver                           |
| QuickGene RNA cultured cell kit S [ for 96 samples ] Reference code RC-S     | ca.10 $\mu$ g / 1 $\times$ 10 <sup>6</sup> cell HL60 cell   |
| QuickGene RNA cultured cell HC kit S [ for 96 samples ] Reference code RC-S2 | 90~150 $\mu$ g / 10cm dish cultured HEK293 cell             |
| QuickGene RNA blood cell kit S [ for 96 samples ] Reference code RB-S        | ca.4.5 $\mu$ g / 1 $\times$ 10 <sup>7</sup> cell Leukocytes |

\*The isolation kits used in QuickGene-Mini80 are also used for QuickGene-810.

## Specifications

### Overview

● Throughput : 1 to 8 samples per run


### Operating conditions

● Supply voltage : AC100~240V  
 ● Power supply frequency : 50 / 60Hz  
 ● Operating conditions : Temperature: 15~30°C  
 Humidity: 30~80%(non-condensing)

### Physical specifications

● Dimensions : 280(W) $\times$ 220(D) $\times$ 180(H)mm  
 ● Weight : Approx. 3kg

\*Research use only



Application Guides Available on our website. ▶ <http://www.kurabo.co.jp/bio/English/>



North American Distributor

[www.autogen.com](http://www.autogen.com)

AutoGen, Inc. 84 October Hill Road Holliston, MA 01746

Tel: 508.429.5965; Fax: 508.429.9765; E-mail: [info@autogen.com](mailto:info@autogen.com)

