


Chapter 3-II-i

Genomic DNA Extraction from Blood of Animal

Genomic DNA Extraction from Bone Marrow Fluid

Protocol



Results


- Electropherogram
No Data
- The yield of genomic DNA
No Data
- Protein contamination : A260/280
No Data
- Chaotropic salt contamination : A260/230
No Data
- Other
No Data

Common protocol is usable for the following

No Data

Genomic DNA Extraction from Buffy Coat

Protocol



*1 Cell number of 3×10^6 were suspended by PBS/200 μ l

Results

■ Electropherogram

No Data

■ The yield of genomic DNA

No Data

■ Protein contamination : A260/280

No Data

■ Chaotropic salt contamination : A260/230

No Data

■ Other


No Data

Common protocol is usable for the following

No Data

Genomic DNA Extraction from Whole Blood of Aves

Protocol



*1 Leave EDB for 30 min at room temperature after adding nuclease-free water and mixing, and use it after perfect dissolution.

Results

■ Electropherogram

No Data

■ The yield of genomic DNA

No Data

■ Protein contamination : A260/280

No Data

■ Chaotropic salt contamination : A260/230

No Data

■ Other


No Data

Common protocol is usable for the following

No Data

Genomic DNA Extraction from Whole Blood of Human

Protocol



*1 Must follow the steps a, b, and c.

*2 Recommend to use the whole blood collected in EDTA-2Na or EDTA-2K.


*3 Proceed the step C immediately after adding whole blood.

*4 Mix completely by vortexing at the maximum speed. If the mixing is not enough by vortexing, use the tapping, pipetting or inverting.

Results

Electropherogram

1 1 1 1 M 2 2 2 2



M : 1k bp ladder
1 : QuickGene
2 : A company (spin method)

The yield of genomic DNA (Sample: 200 μ l of human whole blood)


	(μ g)	Average	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
QuickGene		5.9	7.2	5.3	5.9	5.5	5.5
A company (spin method)		4.5	6.3	4.4	5.2	3.2	3.6

Protein contamination : A260/280

	Average	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
QuickGene	1.94	1.91	1.94	1.96	1.91	1.96
A company (spin method)	1.84	1.86	1.82	1.80	1.87	1.86

Chaotropic salt contamination : A260/230

	Average	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
QuickGene	1.61	1.76	1.69	1.43	1.76	1.42
A company (spin method)	1.12	1.21	0.89	1.07	1.24	1.21




Hemoglobin contamination : A400

	Average	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
QuickGene	0.036	0.023	0.032	0.070	0.031	0.025
A company (spin method)	0.054	0.076	0.040	0.085	0.026	0.043

Other

• PCR

M 1 2 3




Serial dilution of isolated genomic DNA was used for PCR template to amplify p53 exon6 gene.

PCR amplification was performed successfully by using 0.1ng/μl genomic DNA.

M : 100bp ladder
 1 : Genomic DNA 10ng/μl
 2 : Genomic DNA 1ng/μl
 3 : Genomic DNA 0.1ng/μl

• Pulsed-field electropherogram

M1M2 1 2 3




The use of QuickGene-810 (automatic nucleic-acid isolation system) and QuickGene DNA whole blood kit S enables the isolation of long genomic DNA same as manual method using phenol / chloroform.

M1 : MidRange PFG Marker II
 M2 : *Hind* III digest
 1 : Comparison method using spin column (<~70kb)
 2 : Using QuickGene isolation system and reagents (<~140kb)
 3 : Manual method using phenol / chloroform (<~140kb)

• Restriction Enzyme Digestion

1 1 1 1 M 2 2 2 2 M



The eluted genomic DNA sample had been digested with *EcoR* I .
 The success of enzyme digestion is shown by the comparison of lane1 and 2.


M : 1k bp ladder
 1 : Before digestion
 2 : After digestion using *EcoR* I

Common protocol is usable for the following

Canine Whole Blood

Genomic DNA Isolation from Whole Blood of Canine

Protocol



*1 a to c exactly.
Do not add LDB directly after
addition of EDB.

Results

Electropherogram

No Data

The yield of genomic DNA

amount of whole blood	Yield(μ g)
200 μ l	2.52

Protein contamination : A260/280

amount of whole blood	A260/280
200 μ l	1.68

Chaotropic salt contamination : A260/230

amount of whole blood	A260/230
200 μ l	0.61

Other

No Data

Common protocol is usable for the following

Human Whole Blood



North American Distributor
AutoGen, Inc.
84 October Hill Road
Holliston, MA 01746 USA

tel: 508.429.5965
fax: 508.429.9765
email: info@autogen.com
web: autogen.com