

AutoGen Introduces the NeXTRACT 48

AutoGen introduces a new addition to their instrumentation portfolio, the NeXTRACT 48, a fully automated workflow for extracting viral DNA and RNA.

Holliston, MA, December 15, 2020: AutoGen, Inc. is pleased to announce the introduction of their newest fully automated workflow, the NeXTRACT 48. The NeXTRACT 48 is designed for the rapid extraction of viral DNA and/or RNA from a variety of samples including VTM, UVT, and sputum. This workflow processes 48 samples in 20 minutes and can easily process 800 samples or more per day, providing quick, consistent, and reproducible results.

Compatible with most COVID-19 real-time PCR diagnostic kits, the NeXTRACT 48 is compact, inexpensive, and easily fits into a BSL2 environment. The NeXTRACT 48 is user friendly and produces archival quality DNA or RNA with LIMS compatibility. With its small footprint, ease of use, and ability to process 800 samples daily, the NeXTRACT 48 fulfills a need with just one instrument.

“This instrument is a perfect addition for labs looking to supplement COVID-19 viral RNA extractions or to add a back-up platform with low capital investment and operating cost,” said Michael Messier, CEO of AutoGen. “The NeXTRACT 48 provides a solution that customers have been asking for and it is reasonably priced and reliable. Of course, all of this will come with the same unwavering support and industry-leading service that we’ve been providing for over 20 years.”

ABOUT AUTOGEN

AutoGen, a leader in the life sciences marketplace, is an organization that works to understand a lab’s full workflow while identifying areas to improve efficiency and reduce costs within their nucleic acid extraction processes. AutoGen’s workflows provide solutions that are the best fit for their customers’ budgets. Our customers include biorepositories, contract research organizations, academic research laboratories, pharmaceutical companies, clinical diagnostic laboratories, and government institutions all over the world. AutoGen’s mission is to provide quality instrumentation and chemistries, as well as dedicated technical support - all with a level of post-sale service that is truly unmatched.

For more information about AutoGen, visit www.autogen.com. If you would like to learn more about the NeXTRACT 48, contact Rob Osborn at 508-395-8161 or rosborn@autogen.com