

FLEX STAR

AutoGen is the sole provider of precipitation workflows in the marketplace. With 17 years of stability data, the quality of DNA is unmatched.

The FlexSTAR workflow is the only fully automated and cost effective, precipitation-based extraction workflow in the marketplace. Backed by the proven, Qiagen FlexiGene chemistry, no other workflow provides higher quality or higher molecular weight DNA, free of RNA contamination. Upon completion, DNA is ready to use and perfect for applications like long-term storage or straight to any downstream application.



KEY BENEFITS:

- + Isolates DNA from whole blood, buffy coats, tissue, cells, buccal swabs, saliva, and more
- + Load up to 30, 5ml samples or 15, 10ml samples per run.
- + Includes our newly released Oragene® and Isohelix® Saliva protocols allowing flexibility and to process up to 30, 4ml saliva tubes at a time.
- + Includes protocols for buffy coats derived from 10ml (30 at a time) or 20ml (15 at a time)
- + Offers sample tracking integrable with any LIMS system.

SAMPLE TYPE/PROTOCOLS

- + Whole Blood
- + Cultured Cells
- + PBMCs
- + Oragene® Saliva
- + Isohelix® Saiva
- + Fresh/frozen buffy coats
- + Buccal swabs
- + Human/animal tissue



ABOUT AUTOGEN

AutoGen, a leader in the life sciences marketplace, is an organization that works to understand a lab's full workflow while identifying areas to improve efficiency and reduce costs within their nucleic acid extraction processes. AutoGen's workflows provide solutions that are the best fit for their customers' budgets. Our customers include biorepositories, contract research organizations, academic research laboratories, pharmaceutical companies, clinical diagnostic laboratories, and government institutions all over the world. AutoGen's mission is to provide quality instrumentation and chemistries, as well as dedicated technical support - all with a level of post-sale service that is truly unmatched.

For more information about AutoGen, visit www.autogen.com. If you would like to learn more about the FlexSTAR workflows, contact Rob Osborn at 508-395-8161 or rosborn@autogen.com